

Laser Toner Analysis by Pyrolysis-GC/MS

The toner materials used in laser printers and photocopiers are frequently a combination of organic polymers or copolymers and inorganic materials like metal oxides used as pigments. It may be of interest to analyze these toners alone, for manufacturing and quality control purposes, or to analyze them after they have been printed onto paper, as in questioned document investigations. In either case, pyrolysis-GC or pyrolysis-GC/MS offers a simple and effective technique for unraveling such a complex sample.

When paper (which is essentially cellulose) is pyrolyzed, it produces CO_2 , water, oxygenates like furans and aldehydes, small hydrocarbons and levoglucosan. Figure 1 shows a pyrogram of a standard white paper used with laser printers, without toner. When a piece of the same paper which contains toner is pyrolyzed, the pyrogram includes all the cellulose pyrolysate peaks plus peaks from the organic material used in the toner. In this case, the toner was made using a copolymer of styrene and butyl acrylate, and the monomer peaks are easily seen in Figure 2. In fact, the styrene monomer peak is by far the largest peak, and styrene dimer and trimer may be seen as well.

The first ten minutes of these pyrograms have been expanded for easier comparison in Figure 3. Other toner formulations may include additional monomers, such as other acrylics, which are also readily identified using pyrolysis-GC/MS.




Figure 1




Figure 2

Peak Identification for Figures 2 and 3

1. Styrene monomer
2. Butyl acrylate monomer
3. Levoglucosan (from cellulose)
4. Styrene dimer
5. Styrene trimer




Figure 3. Upper - paper only, lower - paper with toner

Equipment

All samples were pyrolyzed using a CDS Model 2500 Pyrolysis Autosampler interfaced to a Hewlett-Packard 6890 gas chromatograph with a mass selective detector.

Pyrolysis

Interface oven: 300°C
 Ramp: 10°C/ms
 Temperature: 750°C
 Time: 15 seconds
 Clean: 1000°C for 10 seconds

Chromatography

Carrier: He
 Column: HP-5
 30 m x 0.25 mm
 Split: 75:1
 Initial temperature: 40°C for 2 min
 Ramp: 10°C/min
 Final temperature: 290°C for 10 min

FOR MORE INFORMATION
 CONCERNING THIS APPLICATION,
 WE RECOMMEND THE
 FOLLOWING READING:

J. Zimmerman, D. Mooney and M. Kimmett, Preliminary Examination of Machine Copier Toners by Infrared Spectrophotometry and Pyrolysis Gas Chromatography, J. Forensic Sci., 31, 2, (1986) 489.

T. Wampler and E. Levy, Applications of Pyrolysis Gas Chromatography/Mass Spectrometry to Toner Materials from Photocopiers, J. Forensic Sci. 31, 1, (1986) 258.

Additional literature on this and related applications may be obtained by contacting your local CDS Analytical representative, or directly from CDS at the address below.



CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 1000, 2000 and 2500 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 6500 16 position autosampler for complex, multicomponent materials investigation. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto www.cdsanalytical.com.