

### Identifying Additives in a Pyrogram using Deconvolution Software

Pyrolysis-GC/MS is a simple technique which may be used to study the composition of complex polymeric materials. Automobile tire rubber, for example is made with polymers, and also additives such as carbon black, antioxidants and other processing agents. The pyrogram of rubber shown in Figure 1 reveals that the polymeric material was formulated using styrene, butadiene and isoprene. Additional information about the additives used is generally also present, requiring the investigation of peaks not related to the rubber pyrolysis products.

Deconvolution software can facilitate this investigation by comparing specific ions in the mass spectra of the various peaks with a library of known compounds. Searching the total ion chromatogram produces a reconstructed TIC made up of the likely matches for compounds in the method library. The peaks are identified by compound name and retention time, making it simple to confirm the identification in the original chromatogram.

Figure 2 shows a deconvolution search of the rubber in Figure 1 performed using the IFD software from Ion Signature and a method library prepared at CDS. At 19.42 minutes, a large peak is indicated for N-(-1-methylethyl)-N'-phenyl-1,4-benzenediamine (called Antioxidant IP among other names). The IFD software indicates the masses used for matching as well as the actual and expected abundances of these masses.




Figure 1. TIC of rubber pyrolyzed at 750°C.




Figure 2. Search results for rubber, indicating Antioxidant IP at 19.42 minutes.

## Equipment

These samples were analyzed using a CDS Model 5200 Pyroprobe, interfaced to an Agilent gas chromatograph/mass spectrometer.

### Model 5200 Conditions

Valve Oven: 300°C  
Transfer Line: 325°C  
Temperature: 750°  
Time: 10 seconds  
Interface: 300°C  
Interface time: 4 min

### GC Conditions

Carrier: Helium  
Injector: 300°C  
Split: 75:1  
Column: 5% phenyl (30m X 0.25mm)  
Detector: Agilent MSD  
Range: 35 - 550

### GC Program:

Initial: 40°C for 2 minutes  
Ramp: 10°C/min.  
Final: 300°C

## Deconvolution Software

Ion Signature

100-01 Quantitative Deconvolution Software

FOR MORE INFORMATION  
CONCERNING THIS APPLICATION,  
WE RECOMMEND THE  
FOLLOWING READING:

Determination of Polymer Additives  
using Analytical Pyrolysis, K. D.  
Jansson et al., J. Anal. Appl. Pyrolysis  
79, (2007) 353-361.

Additional literature on this and related  
applications may be obtained by con-  
tacting your local CDS Analytical rep-  
resentative, or directly from CDS at the  
address below.

CDS Analytical, Inc. has been a leader in the design and manufacture of laboratory instruments for sample preparation and analysis since 1969. We are dedicated to providing the best possible instruments for both research and routine analysis. Well known in the field of pyrolysis, CDS manufactures the Pyroprobe® 5000, 5150, 5200 and 5250 autosampler for the introduction and analysis of solid materials by GC, MS and FT-IR. CDS offers a complete line of dynamic headspace instruments for the analysis of volatile organic compounds in environmental, pharmaceutical and food applications, including the model 8400 four-position autosampler. CDS also manufactures the Dynatherm line of thermal desorption instruments including the 9000 series for air monitoring and the 9300 TDA. Our customers, their requirements and applications are important to us. To help meet your needs, we offer a wide range of analytical information and the services of our applications laboratory. If you would like additional information, please contact us at the address below, call us at 1 800 541 6593, or log onto [www. cdsanalytical.com](http://www.cdsanalytical.com).