

Quantitation of Epoxy to Hardener Ratios in a Cycloaliphatic Epoxy

Application Note

Epoxies

Although many common epoxy resins are formulated using Bisphenol A, aliphatic, especially cyclyhexyl compounds, are also used. The epoxies shown here were formulated with varying amounts of 3,4-Epoxycyclohexylmethyl 3,4- epoxy cyclohexane carboxylate, shown below:



and hardened using methylhexahydrophthalic anhydride (MHHPA), which has the structure:



The samples used epoxy to hardener ratios of 0.50, 0.66 and 1.00. When pyrolyzed, each sample produced a pyrogram like the one shown in Figure 1, with most of the components eluting as three peaks at about 14 minutes. Figure 2 shows an expanded view of the epoxies, with epoxy to hardener ratios of 1:2 and 1:1.

Ascribing the first peak to the epoxy resin and the third peak to the hardener, peak area ratios were plotted agains the relative amounts of the constituents. This produced a linear graph for the range of constituent ratios, shown in Figure 3.



Figure 1. Pyrogram of epoxy at 750°C.

Author:

T. Wampler



Figure 2. Expanded view of pyrogram.

Area ratio vs component ratio



Figure 3. Graph of peak area ratio vs component ratio.

FOR MORE INFORMATION CONCERNING THIS APPLICATION, WE RECOMMEND THE FOLLOWING READING:

H. Nakagawa and S. Tsuge, Studies on Thermal Degradation of Epoxy Resins by High-resolution Pyrolysis-Gas Chromatography, J. Anal. Appl. Pyrolysis 7 (1987) 113.