

Author:

Karen Sam

#245

# **Application Note**

Polymer

## Abstract

This application note investigates the thermal decomposition of Polycarbonate with the presence of steam.

## Introduction

The effects of steam pyrolysis has been studied for both fuel and recycling applications. For example, steam pyrolysis of biomass as an approach for bio-oil production and upgrading has been examined, showing that steam affects the yields and chemical nature of gas, liquid, and char of biomass<sup>1</sup>. On the recycling side, pyrolysis of polycarbonate in the presence of water has been studied, where pyrolysis of polycarbonate can yield some valuable monomer material, but often with too many by-products. On the other hand, hydrolysis of polycarbonate has been shown to increase monomer material<sup>2</sup>. In this application, pyrolysis of Lexan (poly(bisophenol A carbonate)) was studied under different inert and reactive pyrolysis conditions, including temperature ramps, and the presence of steam as a reactant gas.

#### **Experiment Setup**

The sample was first added into a DISC (Drop-In-Sample Chamber) tube and then analyzed using Evolved Gas Analysis (EGA) as a screening step. Using the information from this step, single step pyrolysis was performed. After which, pyrolysis with a 70°C ramp rate was performed both with and without steam as a reactant gas on a Pyroprobe 6200 with Steam Option.

EGA Pyroprobe 6200 DISC: Initial: Ramp Rate: Final: Interface:	100°C 100°C per min 800°C 300°C	GC-MS Column: Carrier: Injector: Oven: Ion Source:	Fused silica (1m x 0.10mm) Helium 1.25mL/min 80:1 split 360°C 300°C 230°C	
Transfer Line:	300°C	Mass Range:	35-600amu	
Valve Oven:	300°C	mass hange.	00-000amu	
Pyrolysis				
Pyroprobe 6200	with Steam	GC-MS		
DISC:	800°C 30s	Column: Carrier:	5% phenyl (30m x 0.25mm) Helium 1.25mL/min	
Interface:	300°C		80:1 split	
Transfer Line:	300°C	Injector:	360°C	
Valve Oven: 300°C		Oven:	40°C for 2 minutes 10°C/min to 320°C(15min)	
		lon Source:	230°C	

Mass Range:

35-600amu



Ramp Rate Pyrolysis with Steam

Pyroprobe 6200 DISC:	with Steam	GC-MS Column:5% phenyl (30m)	
Initial:	300°C 2min	Carrier:	Helium 1.00mL/min
Ramp:	70°C/minute		50:1 split
Final:	800°C 20s	Injector	: 300°C
		Oven:	50°C for 2 min
Purge Gas:	He		15°C/min to 315°C
Flow Rate:	25mL/min	Source:	230°C
		Range:	35-600 amu
Reactant Gas:	Steam		
WaterFlow Rate	: 0.02mL/min		
Trap Contents:	Tenax TA		
Trap Rest:	40°C		
Trap Final:	300°C		
<b>T</b> DI I	2222		
Top Block:	200°C		
Interface:	300°C		
Transfor Line	20000		
Transfer Line:	300°C		
Valve Oven:	325°C		

#### **Results and Discussion**

When subject to a heating rate of 100°C per minute from 100°C to 1000°C, thermal decomposition of polycarbonate begins at 500°C, and is finished by 800°C. So a single step pyrolysis at 800°C was used to initially study the polycarbonate.






Figure 1. Evolved Gas Analysis of Lexan poly(bisphenol A carbonate) from 100°C to 1000°C at 100°C per minute.

At 800°C, Bisphenol A, Lexan's monomer, is a prominent pyrolysis product, along with phenol, cresol, and p-Isopropyl phenol.

Lexan was then pyrolyzed at a 70°C pyrolysis ramp from 300°C to 800°C to a sorbent trap, both with and without steam as a reactant gas. A slow ramp was chosen to allow time for steam to interact with the polymer.



Peak	Identification	Peak	Identification
1	Phenol	6	p-lsopropenylphenol
2	p-Cresol	7	Diphenyl carbonate
3	Phenol, 4-ethyl-	8	Phenol, 4-(1-methyl-1-phenylethyl)-
4	p-Cumenol	9	Bisphenol A
5	Phenol, 4-propyl-		·



When Lexan (poly(bisphenol A carbonate)) is pyrolyzed with a 70°C ramp rate, most of the the same breakdown products are seen but with different relative intensities (Figure 3). When steam is added as a reactant, preliminary results show that bisphenol A is the main product, and much fewer of the other products remain. It is proposed that the main mechanism for production of bisphenol A is the presence of steam allowing for hydrolysis of the carbonate bonds within the polycarbonate polymer.

# Conclusion

A new Steam Option expands the analysis capacity of the Pyroprobe 6200 from CDS Analytical. The reactions with, and the effects of steam on pyrolysis of different materials can be studied.

## References

1. E. Kantarelis, W. Yang, W. Blasiak, Production of Liquid Feedstock from Biomass via Steam Pyrolysis in a Fluidized Bed Reactor., Energy and Fuels, 2013, 27, 4748-4759.

2. G. Bozzano, M. Dente, R. Del Roso, "Poly(bisophenol A carbonate) Recycling: High Pressure Hydrolysis Can be a Convenient Way",Intertech Open, March 2012.